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Wireless Smart Cart [ME-1240 (Red) / ME-1241 (Blue)]


Introduction

The Wireless Smart Cart is a combination wireless and USB device that connects to a computer or tablet through Bluetooth, or to a computer or USB charger through an included micro USB cable. The Smart Cart has a durable ABS body, a three-position plunger, and nearly frictionless wheels. It includes built-in sensors which measure acceleration in three dimensions, position, velocity, and force. The Smart Cart also includes a gyroscope-type sensing element which can measure rotational motion. The Smart Cart can make its measurements on or off a dynamics track and transmit them wirelessly. For easy identification, the Smart Cart is available in red or blue.

The accessory tray on the top of the Smart Cart can hold any accessory designed for a PASCO cart except the Spring Cart Launcher (ME-6843). The Smart Cart has the same dimensions as other PASCO carts and can be used with any PASCO track. The Smart Cart can be stacked with any PASCO cart and has a mass of approximately 0.250 kilograms (250 grams). It has attachment points on the top and bottom of both ends, as well as Velcro® tabs for inelastic collisions on one end. The magnetic bumper, bumper attachment, and hook attachment fit into the threaded hole in the end of the Smart Cart opposite to the tabs.

The built-in sensors measure force in the range between -100 newtons (N) and +100 N, acceleration in the range between -16 g and +16 g (where $g = \pm 9.8 \text{ m/s}^2$), and angular rotation rate up to ± 145 degrees per second (dps). The cart measures force as a push or pull along its x-axis and measures acceleration in three dimensions (x, y, and z). PASCO Capstone or SPARKvue can also show the resultant acceleration. The encoder wheels on the Smart Cart measure velocity to a maximum of 3.0 meters per second; the wheels have a position resolution of 0.2 mm. The positive direction of motion is the same as the x-axis in the printed graphic representing the acceleration sensor's position.


The Smart Cart is designed to optimize battery usage time and will turn off after about five minutes if not connected to the software. Since each Smart Cart has a unique device ID number, more than one can be connected to a computer or tablet at once.

- 1 Three-position plunger
- 2 Acceleration and gyro sensor position
- 3 Attachment point
- 4 Hook attachment
- 5 Attachment point


Equipment

Included equipment:

- 1 Wireless Smart Cart (ME-1240 or ME-1241)
- 2 Magnetic bumper
- 3 Bumper attachment
- 4 Hook attachment
- 5 Micro USB cable (not pictured)

Features:

- 1 Velcro® tab
Allows two carts to attach to each other in an inelastic collision.
- 2 Battery Status LED
Indicates the charge level and status of the cart's battery.

Battery LED	Status
Red blink	Low power
Yellow ON	Charging
Green ON	Fully charged

Wireless Smart Cart | ME-1240 (Red) / ME-1241 (Blue)

3 Bluetooth Status LED

Indicates the status of the cart's Bluetooth connection.

Bluetooth LED	Status
Red blink	Ready to pair
Green blink	Connected

NOTE: The Bluetooth LED will not light up if the sensor is connected to the software via micro USB cable.

4 Device ID

Use this to identify the cart when connecting via Bluetooth.

5 Accelerometer position

Indicates the position of the acceleration and gyro sensor within the cart, as well as the direction of the three dimensions of acceleration (the x-axis, y-axis, and z-axis).

6 Plunger release

Press to immediately release the plunger to its maximum extension.

7 Threaded holes (M5)

Use to secure accessories in place on top of the cart.

8 Accessory tray

Use to hold accessories mounted on the cart.

9 Attachment point

Use these points to tie a string to the cart for various experiments.

10 Micro USB port

Use with the included micro USB cable to connect the cart to a USB charger. The port and cable can also be used to directly connect the cart to a computer without the use of Bluetooth. This connection method is not supported by iOS.

11 Accessory port

Use to connect various accessories to the Smart Cart, enable data collection from them, and control them.

12 ON/OFF button

Press to turn the cart on. Press and briefly hold to turn the cart off.

13 Built-in wheel encoder

Used to automatically record position, velocity, and acceleration of the wheels during data collection.

Required equipment:

- Data collection software (PASCO Capstone or SPARKvue)

Initial step: Charge the battery

Charge the battery by connecting the USB port to any standard USB charger. The Smart Cart Charging Garage (ME-1243) allows up to five Smart Carts to be charged simultaneously, while also serving as a storage option for the cart and bumpers. The battery status light is solid yellow while charging. When fully charged, the light changes to solid green.

Get the software

You can use the Smart Cart with SPARKvue or PASCO Capstone software. If you're not sure which to use, visit [pasco.com/products/guides/software-comparison](https://www.pasco.com/products/guides/software-comparison).

SPARKvue is available as a free app for Chromebook, iOS, and Android devices. We offer a free trial of SPARKvue and Capstone for Windows and Mac. To get the software, go to [pasco.com/downloads](https://www.pasco.com/downloads) or search for **SPARKvue** in your device's app store.

If you have installed the software previously, check that you have the latest update:

SPARKvue

Go to Main Menu  > **Check for Updates**

PASCO Capstone

Go to **Help** > **Check for Updates**.

The Wireless Smart Cart can also be used with the free PASCO program MatchGraph!, which can be downloaded at [pasco.com/matchgraph](https://www.pasco.com/matchgraph).

Check for a firmware update**SPARKvue**

1. Press the power button until the lights turn on.
2. Open SPARKvue.
3. Select **Sensor Data** on the Welcome Screen.

4. Select the sensor that matches your sensor's device ID. A notification appears if a firmware update is available. Click **Yes** to update the firmware.
5. Close SPARKvue.

PASCO Capstone

1. Press and hold the power button until the lights turn on.
2. Open PASCO Capstone.
3. Click **Hardware Setup**.

4. Select the sensor that matches your sensor's device ID. A notification appears if a firmware update is available. Click **Yes** to update the firmware.
5. Close Capstone.

Setting up the software

Collecting data from the Smart Cart requires PASCO Capstone, SPARKvue, or MatchGraph! software. The Smart Cart can be connected to any of these programs using either Bluetooth or a USB connection.

SPARKvue

To connect to SPARKvue using Bluetooth:

1. Turn on the Smart Cart. Check to make sure that the Bluetooth LED is blinking red.
2. Start SPARKvue, then select **Sensor Data** from the main menu.
3. From the list of available wireless devices, select the Smart Cart which matches the device ID (XXX-XXX) printed on your Smart Cart to connect to the cart.

To connect to SPARKvue using the micro USB cable:

1. Turn on the Smart Cart.
2. Start SPARKvue, then select **Sensor Data** from the main menu.
3. Connect the micro USB cable (included) from the micro USB port on the Smart Cart to a USB port on the computer, or into a powered USB hub connected to the computer. The cart will connect automatically.

Upon connecting to the Smart Cart, a list of available Smart Cart sensor measurements will automatically appear in the **Select Measurements for Templates** section.

To collect data with SPARKvue:

1. From the **Sensor Data** menu, after connecting to the Smart Cart, enable the desired measurements by clicking the checkboxes next to those measurements.
2. In the **Templates** section, select **Graph** to enter the Experiment Screen.
3. The Graph template will automatically populate the axes with the measurements you selected on the y-axis (or y-axes if you selected multiple measurements) and time on the x-axis. To change the quantity measured on an axis, click the buttons displaying the current measurement, then select the appropriate quantity from the list.
4. When you are ready, click **Start**  to begin collecting data.

PASCO Capstone

To connect to Capstone using Bluetooth:

1. Turn on the Smart Cart.
2. Start Capstone, then select **Hardware Setup**  from the **Tools** palette.
3. From the list of **Available Wireless Devices**, select the Smart Cart which matches the device ID (XXX-XXX) printed on your Smart Cart to connect to the cart.

To connect to Capstone using the micro USB cable:

1. Turn on the Smart Cart.
2. Start Capstone. If desired, select **Hardware Setup**  from the **Tools** palette to verify the connection status.
3. Connect the micro USB cable (included) from the micro USB port on the Smart Cart to a USB port on the computer, or into a powered USB hub connected to the computer. The cart will connect automatically.

When connecting the Smart Cart to Capstone, the measurement of all sensors will be enabled by default. Any undesired measurements can be disabled through the **Data Summary**  menu in the **Tools** palette.

To collect data from Capstone:

1. Double-click the **Graph** icon in the **Displays** palette to create a new Graph display.
2. On each axis, click **<Select Measurement>**, then choose the appropriate sensor measurement for your experiment from the list.
3. When you are ready, click **Record**  to begin recording data.

MatchGraph!

To connect to MatchGraph! using Bluetooth:

1. Turn on the Smart Cart. Check to make sure that the Bluetooth LED is blinking red.
2. Start MatchGraph!. The **Sensor Not Found** window will open automatically.
3. Select **Choose Wireless Interface** from the bottom right corner of the **Sensor Not Found** window.
4. From the **Select Wireless Interface** list, select the Smart Cart which matches the device ID (XXX-XXX) printed on your Smart Cart to connect to the cart.

To connect to MatchGraph! using the micro USB cable:

1. Turn on the Smart Cart.
2. Connect the micro USB cable (included) from the micro USB port on the Smart Cart to a USB port on the computer, or into a powered USB hub connected to the computer.
3. Start MatchGraph!. The cart will connect automatically.

Note that MatchGraph! is only capable of tracking position and velocity data from the Smart Cart. Data from the cart's other measurements will not be collected.


To collect data from MatchGraph!:

1. Select whether you want to track position or velocity data from the Smart Cart.
2. Choose a profile to be displayed by clicking **View All Profiles**  at the top of the screen and selecting a profile from the list.
3. When you are ready, click **Record**  to begin tracking data.


Wireless Smart Cart | ME-1240 (Red) / ME-1241 (Blue)

Connecting bumper and hook attachments

Screw the bumper or hook attachment into the threaded hole on the front end of the Smart Cart, as shown below.

The magnetic bumper has a metal post that fits in the small hole to the side of the threaded hole on the front end of the Smart Cart, as shown below. The post ensures that the magnetic bumper is mounted in the correct orientation.

Screwing any of these attachments into the threaded hole on the front end of the Smart Cart connects them to the Smart Cart's built-in Force Sensor, allowing measurements of the force on the cart to be collected.

IMPORTANT: If using the Smart Cart for collisions, do NOT use the rubber bumper, as the force of impact on the bumper may damage the internal force sensor. Instead, any of the following components can be used as a bumper for collision experiments:

- Light Spring Bumper (ME-9275)
- Bumper Accessory Set (ME-9884)
- Magnetic Bumper Set (ME-9885A)

Direction of acceleration

The acceleration sensing unit inside the Smart Cart is oriented so that the acceleration axes are aligned with the x-y-z arrows indicating the direction of positive acceleration.

In this orientation, the x-axis direction is 'along' the long axis of the Smart Cart; the y-axis direction is 'perpendicular' relative to the top of the Smart Cart; and the z-axis direction is 'perpendicular' to the long axis of the Smart Cart and parallel relative to the top of the cart.

Zero (tare) the sensors

At the beginning of an experiment, the measurement from the built-in force and acceleration sensors may not be zero when the force or acceleration is actually zero. This is a normal behavior that is corrected when you zero (or *tare*) the sensor using PASCO Capstone or SPARKvue. For information on zeroing a sensor in these programs, see the Capstone or SPARKvue online help and search for "Zero sensor".

Calibration

The Smart Cart is factory calibrated, so calibration is not necessary, especially if you are measuring a change in a measurement rather than absolute values. However, it is possible to calibrate the sensors if desired. For information on calibrating a Smart Cart, see the PASCO Capstone or SPARKvue online help.

NOTE: Zeroing of the acceleration sensor should **only** be performed if you are measuring acceleration on only a single axis, as zeroing the sensor will cause the measurements of acceleration on the other two axes to no longer report accurate values.

Accessories

A number of accessories are available for the Wireless Smart Cart. Links to these products can be found on the product page. In particular, certain accessories are designed specifically to be used with the Wireless Smart Cart, including:

- Smart Fan Accessory (ME-1242)
- Smart Ballistic Cart Accessory (ME-1245)
- Smart Cart Vector Display (ME-1246)
- Smart Cart Motor (ME-1247)

For instructions on mounting these accessories to the Smart Cart, see the manual for the product in question.

Troubleshooting

- If the Wireless Smart Cart loses Bluetooth connection and will not reconnect, try cycling the ON button. Press and briefly **hold** the button until the status LEDs blink in sequence, then release the button. Start the Smart Cart in the usual way.
- If the Smart Cart stops communicating with the computer software or tablet application, try restarting the software or application. If the problem remains, press and hold the ON button for 10 seconds, then release. Start the Smart Cart in the usual way.
- If the above problems do not solve the connection issue, turn Bluetooth off and back on for your computer, then retry.

About the battery

The Smart Cart's battery is partially charged at the factory. If the battery status LED blinks red, use the micro USB cable to connect the sensor to a USB port or a USB charger.

Maximizing battery life

One of the factors that affects battery life is the storage temperature. Therefore, avoid storing the Smart Cart in very cold or very hot environments.

Battery troubleshooting

If the Smart Cart battery will not charge, it may need to be replaced. Contact Technical Support for information about battery replacement.

Software help

The SPARKvue and PASCO Capstone Help provide additional information on how to use this product with the software. You can access the help within the software or online.

SPARKvue

Software: Main Menu  > Help

Online: help.pasco.com/sparkvue

PASCO Capstone

Software: Help > PASCO Capstone Help

Online: help.pasco.com/capstone

Specifications and accessories

Visit the product page at pasco.com/product/me-1240 or pasco.com/product/me-1241 to view the specifications and explore accessories. You can also download experiment files and support documents from the product page.

Experiment files

Download one of several student-ready activities from the PASCO Experiment Library. Experiments include editable student handouts and teacher notes. Visit pasco.com/freelabs/me-1240 or pasco.com/freelabs/me-1241.

Technical support

Need more help? Our knowledgeable and friendly Technical Support staff is ready to answer your questions or walk you through any issues.

-  Chat pasco.com
-  Phone 1-800-772-8700 x1004 (USA)
+1 916 462 8384 (outside USA)
-  Email support@pasco.com

Regulatory information

Limited warranty

For a description of the product warranty, see the Warranty and Returns page at www.pasco.com/legal.

Copyright

This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Trademarks

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.

Product end-of-life disposal

This electronic product is subject to disposal and recycling regulations that vary by country and region.

It is your responsibility to recycle your electronic equipment per your local environmental laws and regulations to ensure that it will be recycled in a manner that protects human health and the environment.

To find out where you can drop off your waste equipment for recycling, please contact your local waste recycle or disposal service, or the place where you purchased the product.

The European Union WEEE (Waste Electronic and Electrical Equipment) symbol on the product or its packaging indicates that this product must not be disposed of in a standard waste container.

CE statement

This device has been tested and found to comply with the essential requirements and other relevant provisions of the applicable EU Directives.

FCC statement

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Battery disposal

Batteries contain chemicals that, if released, may affect the environment and human health.

Batteries should be collected separately for recycling and recycled at a local hazardous material disposal location adhering to your country and local government regulations.

To find out where you can drop off your waste battery for recycling, please contact your local waste disposal service, or the product representative.

The battery used in this product is marked with the European Union symbol for waste batteries to indicate the need for the separate collection and recycling of batteries.

Patents

The Smart Cart (ME-1240 and ME-1241) is covered by U.S. Patent No. 10,481,173.

The Smart Fan Accessory (ME-1242) is covered by U.S. Patent No. 10,482,789.

1.2 Kracht sensor

Product Guide | 012-10237C


High Resolution Force Sensor (PS-2189)

Introduction

The High Resolution Force Sensor connects to a PASPORT interface and records force in the range of -50 N to +50 N at a rate of up to 1000 samples per second. High resolution, dynamic over-sampling, and low drift rate make this sensor well-suited for measuring very small changes in force.

Equipment

Included equipment:

- ❶ High Resolution Force Sensor
- ❷ Bumper attachment
- ❸ Hook attachment
- ❹ Cart thumbscrew (M5 × 45 mm)
Use to secure the sensor to a PASCO cart.
- ❺ Rod clamp thumbscrew (1/4-20 × 0.75 in.)
Use to secure the sensor to a rod, such as the 120 cm Stainless Steel Rod (ME-8741).

Additional items required:

- PASPORT interface
- PASCO Capstone or SPARKvue data collection software

Get the software

You can use the sensor with SPARKvue or PASCO Capstone software. If you're not sure which to use, visit [pasco.com/products/guides/software-comparison](https://www.pasco.com/products/guides/software-comparison).

SPARKvue is available as a free app for Chromebook, iOS, and Android devices. We offer a free trial of SPARKvue and Capstone for Windows and Mac. To get the software, go to [pasco.com/downloads](https://www.pasco.com/downloads) or search for **SPARKvue** in your device's app store.

If you have installed the software previously, check that you have the latest update:

📱 **SPARKvue**

Go to Main Menu > Check for Updates

📁 **PASCO Capstone**


Go to Help > Check for Updates.

Set up the hardware

Sensor mounting


To mount the sensor on a rod:

1. Slide the sensor onto a rod, as shown below.
2. Tighten the thumbscrew to secure the sensor in place.

To mount the sensor on a PASCO cart:

1. Align the hole labeled **Cart** on the sensor with one of the threaded holes in the accessory tray of the PASCO cart.
2. Insert the included cart thumbscrew through the **Cart** hole so that it extends into the threaded hole, as shown below.
3. Tighten the thumbscrew to secure the sensor in place.


PASCO

1

High Resolution Force Sensor | PS-2189

Connect bumper and hook attachments

Screw the bumper or hook into the sensor as shown below.

IMPORTANT: When performing experiments that involve a collision, **do NOT** use the rubber bumper with the Force Sensor, as this may damage the internal components. Instead, any of the following components can be used as a bumper for collision experiments:

- Light Spring Bumper (ME-9275)
- Bumper Accessory Set (ME-9884)
- Magnetic Bumper Set (ME-9885A)

Set up the software **SPARKvue****Connecting the sensor to SPARKvue:**

1. Turn on SPARKvue, then click **Sensor Data**.
2. Turn on the PASPORT interface if needed, then connect your interface to SPARKvue. For more specific details, see the manual for your chosen interface and the SPARKvue online help.
3. Plug the High Resolution Force Sensor into one of the PASPORT ports on the interface. SPARKvue will automatically detect and identify the sensor.

Collecting data:

1. From the **Select Measurements for Templates** column, select the appropriate measurement for your experiment.
2. From the **Templates** column, click **Graph** to enter the experiment screen. The display will automatically plot your selected measurement on the y-axis and time on the x-axis.
3. When you are ready, click **Start**  to begin recording data.

 **PASCO Capstone****Connecting the sensor to Capstone:**

1. Turn on PASCO Capstone, then click **Hardware Setup** in the **Tools** palette.
2. Turn on the PASPORT interface if needed, then connect your interface to Capstone. For more specific details, see the manual for your chosen interface and the Capstone online help.
3. Plug the High Resolution Force Sensor into one of the PASPORT ports on the interface. Capstone will automatically detect and identify the sensor.

Collecting data:

1. Double-click the **Graph** icon in the **Displays** palette to create a Graph display.
2. Click each **<Select Measurement>** box and select an appropriate measurement for your experiment to assign that measurement to the associated axis.
3. When you are ready, click **Record**  to begin recording data.

About the force measurement**Zeroing the sensor**

Press the **ZERO** button on the sensor to adjust the sensor's output to zero. The sensor should always be zeroed before an experiment to ensure accurate values are collected.

Sample rate

By default, the sensor collects 20 samples per second. The sensor can collect data as fast as 1000 samples per second (or 2000 samples per second if using an 850 or 550 Universal Interface) and as slowly as one sample every 24 hours. The sample rate can be changed in the data collection software.

Dynamic over-sampling

Over-sampling reduces noise, produces slower data, and improves the measurement resolution. This effect is especially noticeable when very small force changes are measured. The degree of over-sampling depends on the sample rate. To increase the degree of dynamic over-sampling, lower the sample rate. Maximum over-sampling occurs at sample rates of 20 Hz and slower.

Software help

The SPARKvue and PASCO Capstone Help provide additional information on how to use this product with the software. You can access the help within the software or online.

 **SPARKvue**

Software: Main Menu  > Help

Online: help.pasco.com/sparkvue

 **PASCO Capstone**

Software: Help > PASCO Capstone Help

Online: help.pasco.com/capstone

Specifications and accessories

Visit the product page at pasco.com/product/PS-2189 to view the specifications and explore accessories. You can also download experiment files and support documents from the product page.

Experiment files

Download one of several student-ready activities from the PASCO Experiment Library. Experiments include editable student handouts and teacher notes. Visit pasco.com/freelabs/PS-2189.

Technical support

Need more help? Our knowledgeable and friendly Technical Support staff is ready to answer your questions or walk you through any issues.

-  Chat [pasco.com](https://www.pasco.com)
-  Phone 1-800-772-8700 x1004 (USA)
+1 916 462 8384 (outside USA)
-  Email support@pasco.com

Regulatory information

Limited warranty

For a description of the product warranty, see the Warranty and Returns page at www.pasco.com/legal.

Copyright

This document is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Trademarks

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.

Product end-of-life disposal

This electronic product is subject to disposal and recycling regulations that vary by country and region. It is your responsibility to recycle your electronic equipment per your local environmental laws and regulations to ensure that it will be recycled in a manner that protects human health and the environment. To find out where you can drop off your waste equipment for recycling, please contact your local waste recycle or disposal service, or the place where you purchased the product. The European Union WEEE (Waste Electronic and Electrical Equipment) symbol on the product or its packaging indicates that this product must not be disposed of in a standard waste container.


1.3 Rotatie sensor

PASCO

Reference Guide
013-15949A

Wireless Rotary Motion Sensor

PS-3220

Equipment Included	Equipment Included
Wireless Rotary Motion Sensor	Three-step Pulley
O-ring	Small Thumbscrew
Micro USB Cable (1 meter)	Thumbscrew

Required Item*	Part Number
PASCO Data Collection Software: Capstone or SPARKvue	see www.pasco.com

*See the PASCO catalog or the PASCO web site for more information.

www.pasco.com

For downloadable experiments, go to www.pasco.com and enter PS-3220 in the Search window. Check under Resources.

+1 916 462 8384
support@pasco.com


800-772-8700
www.pasco.com

Equipment used with the sensor	Equipment used with the sensor
Rotational Inertia Accessory Kit (ME-3420) (includes Ring and Disk Set (ME-3419), Pendulum Accessory (ME-8969), and Super Pulley with Clamp (ME-9444B))	
Physical Pendulum Set (ME-9633)	Track String Adapter (ME-6569)
Rotating Platform (ME-8951)	Centripetal Force Pendulum (ME-9821)
Gyroscope Mounting Bracket (ME-8963)	Polarization Analyzer (OS-8533A)
Dynamics Track Mount Accessory (CI-6692)	Three-Axis Gyroscope (ME-8960)
Three-Step Pulley Accessory (CI-6693)	

Introduction

The PASCO Wireless Rotary Motion Sensor is a versatile position and motion measuring device. It measures angles to a resolution of 0.18° , and detects the direction of motion. Markings on the outside of the case indicate which is the default positive direction. The maximum speed is between 20 and 80 revolutions per second.

The sensor comes with a removable three-step pulley (10 millimeter (mm), 29 mm, and 48 mm diameters), and a rubber "O"-ring that fits into the largest diameter step. The three-step pulley can be placed large-diameter down or large-diameter up on the shaft. A tab on the inside of the pulley matches a notch on the outside of the shaft. The pulley has a notch and a small hole in the outer edge of the largest and second largest steps for attaching a string. One side of the sensor has a platform for mounting a Super Pulley with Table Clamp (optional).

Using the Sensor

The sensor is designed to work with PASCO data collection software to measure position, velocity, and acceleration. Use it to study optics, dynamics, centripetal force, or motion of a pendulum.

Data Collection Software

PASCO Capstone

- Mac OS X
- Windows

SPARKvue

- Mac OS X
- Windows
- iOS
- Android
- Chromebook

See the PASCO web site at

www.pasco.com/software

for help in selecting the right PASCO software and to check for the latest versions.

Software Help

See the SPARKvue Help or PASCO Capstone Help for information about collecting, displaying, and analyzing data.

- In SPARKvue, select the HELP button (?) in any screen including the Home Screen.
- In PASCO Capstone, select **PASCO Capstone Help** from the **Help** menu, or press **F1**.

BluetoothSM Compatibility

For more information about wireless compatibility, see the PASCO website at:


www.pasco.com/compatibility

Platform	Bluetooth SMART Compatibility
iOS	iPad 3 and later iPhone 4S and later iPod touch 5 and later
SPARK LX / LXi	All models
Android	Android 4.4 and later
Chromebook	Chrome OS (requires PS-3500 Adapter*)
Mac OS X	Models introduced July 2011 or later*
Windows	Windows 7 and later (requires PS-3500 Adapter*)

See Appendix B for more information about the PS-3500 Adapter and Mac OS X models.

Initial Step: Charge the Battery

- **Connect the Cable:** Use the Micro USB Cable to connect the micro USB port on the side of the sensor to a USB **port** or USB **charger** such as the PASCO PS-2575 USB Single Port Charger. Charging begins automatically. The charger circuit inside the sensor turns itself off when the unit is fully charged. The battery status LED will shine yellow as the battery is charging, and will shine green when the battery is charged. The battery is partially charged at the factory. Initial charging time may be three hours or longer depending on the power source and the condition of the battery.

ON/OFF Information

To turn the sensor on, press the ON button. The status LEDs will blink. To turn the sensor off, press and **hold** the ON button for a moment until the status LEDs stop blinking. The sensor puts itself to sleep after one hour of inactivity if connected, and after several minutes if not connected.

Wireless Rotary Motion Sensor

Set Up the Software

LED Information

The Bluetooth and the Battery Status LEDs operate as follows:

For a wireless Bluetooth connection:

Bluetooth LED	Status
Red blink	Ready to pair
Green blink	Connected
Yellow blink	Logging*

Battery LED	Status
Red blink	Low power

For a micro USB cable connection to a USB port:

Bluetooth LED	Status
OFF	--
OFF	--
Yellow blink	Logging*

Battery LED	Status
Yellow ON	Charging
Green ON	Charged

For a micro USB cable connection to a USB charger:

Bluetooth LED	Status
Red blink	Ready to pair
Green blink	Connected
Yellow blink	Logging*

Battery LED	Status
Yellow ON	Charging
Green ON	Charged

***Logging:** PASCO wireless sensors can either stream live data to a compatible device or log data (save it to the sensor's memory). The data can then be uploaded to the device for display and analysis at a later time. Logging capability supports long-term or remote data collection while not connected to the device.

Note: The latest versions of SPARKvue and PASCO Capstone support logging. Check the PASCO Web page at:

www.pasco.com/software

for the latest software version.

Set Up the Software**SPARKvue****Connecting the Sensor to a Tablet or a Computer via Bluetooth**

- For SPARKvue, select the Bluetooth icon. In the **Wireless Devices** list. The sensors are ordered by proximity to the device. Select the correct address that matches the Device ID XXX-XXX number found on the sensor. Select **Done**.

Connecting the Sensor to a Computer with the Micro USB Cable

- Connect the micro end of the included Micro USB Cable into the micro USB port on the side of the sensor. Connect the other end of the Micro USB Cable to a USB port on the computer, or into a powered USB hub connected to the computer.
- In the SPARKvue Home Screen, select a measurement from the list under the sensor's name. A graph of the measurement versus time opens.

PS-3220

Mounting the Wireless Rotary Motion Sensor

Collecting Data

- Select the Start button to begin collecting data.

PASCO Capstone**Connecting the Sensor to a Tablet or a Computer via Bluetooth**

- For PASCO Capstone, select **Hardware Setup** in the Tools palette. In **Hardware Setup** the sensors are ordered by proximity to the device. Select the address that matches the Device ID XXX-XXX number on the sensor.

Select a display in the main window or from the **Display** palette. In the display, use the **<Select Measurement>** menu to pick a measurement to be shown.

Connecting the Sensor to a Computer with the Micro USB Cable

- Connect the micro end of the included Micro USB Cable into the micro USB port on the end of the sensor. Connect the other end of the Micro USB Cable to a USB port on the computer, or into a powered USB hub connected to the computer.
- In PASCO Capstone, select a display in the main window or from the **Displays** palette. In the display, use the **<Select Measurement>** menus to pick the measurement to be shown.

Collecting Data

- Select **Record** to begin recording data.

Troubleshooting the Sensor


- If the sensor loses Bluetooth connection and will not reconnect, try cycling the ON button. Press and briefly **hold** the button until the status LEDs blink in sequence, and then release the button. Start the sensor in the usual way.
- If the sensor stops communicating with the computer software or tablet application, try restarting the software or application. If the problem remains, press and **hold** the ON button for 10 seconds and then release. Start the sensor in the usual way.
- Turn Bluetooth off and then back on. Retry.

Mounting the Wireless Rotary Motion Sensor

The sensor case has two support rod holes that fit rods up to 12.7 mm in diameter, such as the ME-8736 45 cm Stainless Steel Rod, and the case can be put on the support rod using either of the support rod holes.

It is possible to mount the Wireless Rotary Motion Sensor horizontally on a support rod with the Three-step Pulley facing up or facing sideways. You can mount the sensor vertically with the pulley facing forward.

When mounted on a track as shown, a Rotary Motion Sensor could be used to measure the motion of a PASCO Cart as it is pulled by a string suspended over the Three-step Pulley of the sensor and attached to a hanging mass.

Wireless Rotary Motion Sensor Mounting the Wireless Rotary Motion Sensor


Attaching Accessories to the Wireless Rotary Motion Sensor

Using the Disk from the Disk and Ring Set (ME-3419)

For rotational inertia experiments, mount the Wireless Rotary Motion Sensor with the Disk horizontal. The underside of the Disk has a square shaped indent that fits the square-shaped top of the 3-Step Pulley. Mount a Super Pulley with Clamp (ME-9448B) on the platform at the side of the sensor. To provide a known torque, wrap a string around one of the steps of the 3-step Pulley. Arrange the end of the string over the Super Pulley and attach a hanging mass.

Adjust the Super Pulley on the platform so that the string is tangent to the step on the 3-step Pulley on the sensor.

Perform a conservation of angular momentum experiment by dropping a second Disk onto the first Disk as it is rotating. Alternately, attach one of the Ring Alignment Devices to the first Disk and then drop the Ring onto the first disk. (See the suggested experiments.)


PS-3220

Mounting the Wireless Rotary Motion Sensor

Attaching the Pendulum Accessory Rod (part of ME-8969 Pendulum Accessory) to the Wireless Rotary Motion Sensor

To mount the rod of the Pendulum Accessory to the Wireless Rotary Motion Sensor, orient the 3-step Pulley so that the large diameter step has two pair of rod guides opposite each other on the top edge. Align the rod with the rod guides and use the captive screw in the center of the rod to attach the rod and pulley onto the sensor's shaft.

Point Mass Setup Attach the center of the rod to the 3-step Pulley and shaft and mount the cylindrical masses at the ends of the rod to investigate the rotational inertia (moment of inertia) of point masses.


Pendulum Setup Attach the end of the rod to the 3-step Pulley and shaft. Mount a cylindrical mass on the rod to use the rod as a pendulum. Investigate the period of oscillation of the pendulum when the amount of mass or the position of the mass is changed. Investigate the period of oscillation as the amplitude of the swing is changed.

Mounting the Wireless Rotary Motion Sensor on a PASCO Track

The sensor can also be mounted on the short rod that is part of the Dynamics Track Mount Accessory (CI-6692).


Slide the square nut of the Dynamics Track Mount Accessory into the T-slot on the side of the track. Adjust the position of the mounting rod on the Dynamics Track Mount Accessory.

When mounted on the track as shown, the Wireless Rotary Motion Sensor could be used to measure the motion of a PASCO Cart as it is pulled by a string suspended over the Three-step Pulley of sensor and attached to a mass hanger.

Wireless Rotary Motion Sensor

Experiments


See the PASCO web site at www.pasco.com for more information.

Mounting the Wireless Rotary Motion Sensor to a Gyroscope

Use the Gyroscope Mounting Bracket (ME-8963, available separately) to mount the Wireless Rotary Motion Sensor to the Three-Axis Gyroscope (ME-8960). Remove the gyroscope assembly from the vertical shaft of the large "A" base. Mount the Gyroscope Mounting Bracket onto the vertical shaft. Remove the three-step pulley from the Wireless Rotary Motion Sensor, and mount the sensor on the Gyroscope Mounting Bracket with the thumbscrews included with the bracket. Mount the slotted guide arm onto the shaft of the Wireless Rotary Motion Sensor. Replace the gyroscope assembly on the vertical shaft.

See the PASCO web site at www.pasco.com for more information.

**Experiments**

Experiments for the Wireless Rotary Motion Sensor in electronic format are available to download from the PASCO web site.

www.pasco.com

Go to the web site, enter PS-3220 in the Search window, and check under Resources.

Three suggested experiments are:

- Rotational Inertia of a Point Mass
- Rotational Inertia of Disk and Ring
- Conservation of Angular Momentum.

Specifications

Item:	Value
Three-step Pulley	10, 29 and 48 mm diameter
Sensor Dimensions	9 cm by 6.5 cm by 4 cm, 6.35 mm diameter shaft
Resolution	$\pm 0.09^\circ$ or 0.0078 mm 0.02 mm (linear) and 0.09° (angular) at 2,000 points per revolution

PS-3220

Experiments

Item:	Value
Rotational Resolution	0.00157 radian
Maximum Rotation Rate	30 rotations per second
Optical Encoder	Bidirectional, indicates direction of motion, 4,000 divisions/rotation

Technical Support

For assistance with any PASCO product, contact PASCO at:

Address: PASCO scientific
10101 Foothills Blvd.
Roseville, CA 95747-7100

Phone: +1 916-462-8384 (worldwide)
877-373-0300 (U.S.)

Web: www.pasco.com

Email: support@pasco.com

Check the PASCO website for the latest version of the instruction manual.

www.pasco.com/manuals

Limited Warranty For a description of the product warranty, see the PASCO catalog. **Copyright** The PASCO scientific *Instruction Manual* is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited. Rev: 11/18. **Trademarks** PASCO, PASCO Capstone, and SPARKvue are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. For more information visit

www.pasco.com/legal.

Product End of Life Disposal Instructions:

This electronic product is subject to disposal and recycling regulations that vary by country and region. It is your responsibility to recycle your electronic equipment per your local environmental laws and regulations to ensure that it will be recycled in a manner that protects human health and the environment. To find out where you can drop off your waste equipment for recycling, please contact your local waste recycle/disposal service, or the place where you purchased the product.

The European Union WEEE (Waste Electronic and Electrical Equipment) symbol (to the right) and on the product or its packaging indicates that this product **must not** be disposed of in a standard waste container.

Appendix A: BluetoothSM Compatibility

Check the PASCO Web page at

www.pasco.com/compatibility

for the latest information on Bluetooth SMART compatibility.

Platform	Bluetooth SMART Compatibility
iOS	iPad 3 and later iPhone 4S and later iPod touch 5 and later
SPARK Element	All models
Android	Android 4.3 and later
Chromebook	Chrome OS (requires PS-3500 Adapter*)
Mac OS X ¹	Models introduced July 2011 or later
Windows 7 and 8	Requires PS-3500 Adapter*
Windows 10	Bluetooth SMART compatible

PS-3220

Appendix A: Compatibility

*The PS-3500 USB Bluetooth 4.0 Adapter, when connected to a USB port, allows up to three Bluetooth SMART devices, such as this PASCO wireless device, to connect to Windows computers, Chromebooks, and older Macintosh computers.

Note: The PS-3500 USB Bluetooth 4.0 Adapter is the only adapter we can currently recommend. Many other Bluetooth 4.0 adapters are available but this adapter has a specific design that enables in-app pairing of Bluetooth SMART sensors.

¹To check the Mac computer's Bluetooth compatibility, do the following:

- Click the  (Apple) Menu.
- Select *About This Mac*
- Click the *More Info...* button.
- Click the *System Report...* button.
- Select *Bluetooth* from the sidebar on the left, underneath *Hardware*.
- Scan down the list of information until you find "LMP Version".
- If your Mac is equipped with Bluetooth SMART, the LMP Version will show **0x6**. (Anything lower than **0x6** means an older version of Bluetooth. Your device will need the PS-3500 USB Bluetooth 4.0 Adapter.)

¹The Mac Mini and MacBook Air were updated with Bluetooth SMART support in 2011. The MacBook Pro was updated in 2012. The Mac Pro that debuted in December 2013 has Bluetooth SMART support.

Exception: Before you upgrade to El Capitan (Mac OS X 10.11.x), if you have a Macintosh with LMP version "0x4" that requires the PS-3500 USB Bluetooth 4.0 Adapter, please contact PASCO Technical Support for further instructions.

Experiment 1: Rotational Inertia of a Point Mass

Equipment Required*	Equipment Required*
Wireless Rotary Motion Sensor (PS-3220)	Base and Support Rod (ME-9355)
PASCO Data Collection Software	Mass and Hanger Set (ME-8979)1
Rotational Inertia Accessory Kit (ME-3420)	Triple Beam Balance (SE-8723)
Calipers (SF-8711)	Paper clips (for masses <1g)

*See the PASCO Web site at www.pasco.com for more information

Purpose

The purpose of this experiment is to find the rotational inertia of a point mass experimentally and to verify that this value corresponds to the calculated theoretical value.

Theory

Theoretically, the rotational inertia, I , of a point mass is given by $I = MR^2$, where M is the mass, and R is the distance the mass is from the axis of rotation. Since this experiment uses two masses equidistant from the center of rotation, the total rotational inertia will be

$$I_{total} = M_{total}R^2$$

where $M_{total} = M_1 + M_2$, the total mass of both point masses.

To find the rotational inertia experimentally, a known torque is applied to the object and the resulting angular acceleration is measured. Since $\tau = I\alpha$,

$$I = \tau/\alpha$$

where α is the angular acceleration, which is equal to a/r (a = linear acceleration), and τ is the torque caused by the weight hanging from the thread that is wrapped around the 3-step Pulley.

$$\tau = rT$$

where r is the radius of the chosen pulley about which the thread is wound, and T is the tension in the thread when the apparatus is rotating.

Applying Newton's Second Law for the hanging mass, m , gives

$$\Sigma F = mg - T = ma$$

(see Figure 1.1). Solving for the tension in the thread gives:


$$T = m(g - a)$$

After the angular acceleration of the mass (m) is measured, the torque and the linear acceleration can be obtained for the calculation of the rotational inertia.

Wireless Rotary Motion Sensor Experiment 1: Rotational Inertia of a Point

Equipment Setup

1. Attach a mass on each end of the rod (part of the Rotational Inertia (Accessory Kit) equidistant from the rod center. You may choose any radius you wish.
2. Tie one end of a thread to a Mass Hanger and tie the other end to one of the levels of the 3-step Pulley on the Wireless Rotary Motion Sensor (WRMS).
3. Mount the rod and masses to the pulley on the Wireless Rotary Motion Sensor. Please note the orientation of the 3-step Pulley.
4. Mount the WRMS on a support rod. Make sure that the support rod does not interfere with the rotation of the rod and masses. See Figure 1.1.
5. Mount the Super Pulley with Clamp (part of the ME-3420) on the platform of the Wireless Rotary Motion Sensor.
6. Drape the thread over the Super Pulley such that the thread is in the groove of the pulley and the Mass Hanger hangs freely (see Figure 1.1).

Note: The Super Pulley with Clamp must be adjusted at an angle, so that the thread runs in a line tangent to the point where it leaves the 3-step Pulley and straight down the middle of the groove on the clamp-on Super Pulley.

7. Adjust the Super Pulley height so that the thread is level with the 3-step pulley.

Procedure

Part 1: Measurements for the Theoretical Rotational Inertia

1. Weigh the two masses from the ends of the thin rod to find the total mass M_{total} and record the value in Data Table 1.
2. Measure the distance from the axis of rotation to the center of the masses and record this radius in Data Table 1.

Data Table 1: Theoretical Rotational Inertia

Total Mass	
Radius	

Part 2: Measurement for the Experimental Method

A. Finding the Acceleration of the Point Masses and Apparatus

1. In the data collection software, create an experiment to measure the angular velocity (in radians per second) versus time (in seconds) of the point masses and apparatus.
 - In PASCO Capstone, for example, drag the Graph icon from the Displays palette to the workbook. Select "Angular Velocity (rad/s)" for the vertical axis, and "Time (s)" for the horizontal axis.
 - Click the Hardware Setup icon in the Tools palette to open the "Hardware Setup" panel. In the panel, click the properties button (it looks like a gear wheel in the lower right corner).

PS-3220

Experiment 1: Rotational Inertia of a Point Mass

- In the Properties window for Linear Accessory, select the size of the 3-step Pulley you are using. The default setting is "Large Pulley (Groove)". Click OK.
2. Put a 50-g mass on the Mass Hanger and turn the 3-step Pulley to wind up the thread so the hanger is just below the Super Pulley. Hold the 3-step Pulley.
3. Click Record to begin recording data, and release the 3-step Pulley, allowing the hanger to fall.
4. Caution! Click Stop to end data recording BEFORE the hanger reaches the floor or the thread completely unwinds from the 3-step Pulley.
5. In the Graph display, select the region of the data that represents when the Point Masses and Apparatus were accelerating.
6. In the display, select "Linear" from the curve fit menu.

The slope, m , of the linear fit represents the angular acceleration (a) for the Point Mass and Apparatus

7. Record the value of the slope, m , as the angular acceleration in Data Table 2.
8. Using calipers, measure the diameter of the pulley about which the thread is wrapped and calculate the radius. Record the radius in Data Table 2.

In the previous procedure, the apparatus is rotating and contributing to the total rotational inertia. The next step is to find the rotational inertia of the apparatus by itself so that this rotational inertia can be subtracted from the total.

Wireless Rotary Motion Sensor Experiment 1: Rotational Inertia of a Point

B. Finding the Acceleration of the Apparatus Alone

1. Take the point masses off the ends of the rod.
2. Repeat the procedure from Part A for finding the angular acceleration of the apparatus alone.
 - You may need to decrease the amount of hanging mass so that the apparatus does not accelerate too fast for smooth data collection.
 - Remember that the value of the slope, m , is the angular acceleration.
3. Record the data in Data Table 2.

Data Table 2: Experimental Rotational Inertia Data

	Point Masses and Apparatus	Apparatus Alone
Hanging Mass		
Slope, m		
Radius		

Calculations

1. Calculate the experimental value of the rotational inertia of the point masses and apparatus together and record the calculation in Data Table 3.
2. Calculate the experimental value of the rotational inertia of the apparatus alone and record the calculation in Data Table 3.
3. Subtract the rotational inertia of the apparatus from the total rotational inertia of the point masses and apparatus together. Record this in Data Table 3 as the rotational inertia of the point masses alone.
4. Calculate the theoretical value of the rotational inertia of the point masses and record the calculation in Data Table 3.
5. Calculate the percent difference to compare the experimental value to the theoretical value, and record the percent difference in Data Table 3.

Data Table 3: Results

Component	Rotational Inertia
Point Masses and Apparatus Combined	
Apparatus Alone	
Point Masses (experimental value)	
Point Masses (theoretical value)	
Percent Difference	

Experiment 2: Rotational Inertia of Disk and Ring

Equipment Required*

Wireless Rotary Motion Sensor (PS-3220)
 PASCO Data Collection Software
 Rotational Inertia Accessory Kit (ME-3420)
 Calipers (SF-8711)

Equipment Required*

Base and Support Rod (ME-9355)
 Mass and Hanger Set (ME-8979)1
 Triple Beam Balance (SE-8723)
 Paper clips (for masses <1g)

*See the PASCO Web site at www.pasco.com for more information

Purpose

The purpose of this experiment is to experimentally find the rotational inertia of a ring and a disk and verify that these values correspond to the calculated theoretical values.

Theory

Theoretically, the rotational inertia, I , of a ring about its center of mass is given by:

$$I = \frac{1}{2}M(R_1^2 + R_2^2)$$

where M is the mass of the ring, R_1 is the inner radius of the ring, and R_2 is the outer radius of the ring. See Figure 2.1.

The rotational inertia of a disk about its center of mass is given by:

$$I = \frac{1}{2}MR^2$$

where M is the mass of the disk and R is the radius of the disk. See Figure 2.2. To find the rotational inertia experimentally, a known torque is applied to the object and the resulting angular acceleration is measured. Since $\tau = I\alpha$,

$$I = \frac{\tau}{\alpha}$$

where α is the angular acceleration, which is equal to a/r (a = acceleration), and τ is the torque caused by the weight hanging from the thread that is wrapped about the 3-step Pulley on the Rotary Motion Sensor. The torque is given by:

$$\tau = rT$$

where r is the radius of the pulley step about which the thread is wound, and T is the tension in the thread when the apparatus is rotating.

Applying Newton's Second Law for the hanging mass, m , gives:

$$\Sigma F = mg - T = ma$$

See Figure 2.3. Solving for the tension in the thread gives:

$$T = m(g - a)$$

Once the angular acceleration is measured, the radius and the linear acceleration, a , can be obtained for the calculation of the torque.


Figure 2.1: Ring about center of mass


Figure 2.2: Disk about center of mass

Wireless Rotary Motion Sensor Experiment 2: Rotational Inertia of Disk and

Procedure

Measurements for the Theoretical Rotational Inertia


1. Weigh the ring and the disk to find their masses and record these masses in Data Table 1.
2. Measure the inside and outside diameters of the ring and calculate the radii, R_1 and R_2 . Record in Data Table 1.
3. Measure the diameter of the disk and calculate the radius, R , and record into Data Table 1.

Data Table 1: Theoretical Rotational Inertia.

Mass of ring	
Mass of disk	
Inner radius of ring	
Outer radius of ring	
Radius of disk	

Setup

1. Mount the Wireless Rotary Motion Sensor (WRMS) on a support rod.
2. Mount the Super Pulley with Clamp on the side of the WRMS as shown in Figure 2.3.
3. Tie one end of a thread to a Mass Hanger and the other end of the thread to one of the levels of the 3-step Pulley on the RMS.
4. Drape the thread over the Super Pulley such that the thread is in the groove of the Super Pulley and the Mass Hanger hangs freely.
5. Adjust the Super Pulley with Clamp to an angle so that the thread runs in a line tangent to the point where it leaves the 3-step Pulley and is straight down the middle of the groove on the Super Pulley.
6. Remove the thumbscrew from the 3-Step Pulley. Place the disk directly on the square-shaped top of the 3-step Pulley as shown in Figure 2.3.
7. Place one of the disk alignment guides on the disk. Replace the thumbscrew to the 3-Step Pulley so that the disk alignment guide is firmly held in place.
8. Place the thin ring on the disk alignment guide on top of the disk as shown in figure 2.3..


PS-3220

Experiment 2: Rotational Inertia of Disk and Ring

Procedure**Measurements for the Experimental Method****A. Finding the Acceleration of the Ring and Disk**

1. In the data collection software, create an experiment to measure the angular velocity (in radians per second) versus time (in second) of the ring and disk.
 - In PASCO Capstone, for example, drag the Graph icon from the Displays palette to the workbook. Select "Angular Velocity (rad/s)" for the vertical axis, and "Time (s)" for the horizontal axis.
 - Click the Hardware Setup icon in the Tools palette to open the "Hardware Setup" panel. In the panel, click the properties button (it looks like a gear wheel in the lower right corner).
 - In the Properties window for Linear Accessory, select the size of the 3-step Pulley you are using. The default setting is "Large Pulley (Groove)". Click OK.
2. Put a 50-g mass on the Mass Hanger and turn the 3-step Pulley to wind up the thread so the hanger is just below the Super Pulley. Hold the 3-step Pulley.
3. Click Record to begin recording data, and release the 3-step Pulley, allowing the hanger to fall.
4. Caution! Click Stop to end data recording BEFORE the hanger reaches the floor or the thread completely unwinds from the 3-step Pulley.
5. In the Graph display, select the region of the data that represents when the ring and disk were accelerating.
6. In the Graph display, select "Linear" from the curve fit menu.

The slope, m , of the linear fit represents the angular acceleration (α) for the Point Mass and Apparatus

7. Record the value of the slope, m , as the angular acceleration in Data Table 2.

Data Table 2: Experimental Rotational Inertia Data

	Ring and Disk Combined	Disk Alone
Hanging Mass		
Slope, m		
Radius		

C. Finding the Acceleration of the Disk Alone

1. In "Finding the Acceleration of Ring and Disk," both the disk and the ring are rotating; therefore, it is necessary to determine the acceleration and the rotational inertia of the disk by itself so this rotational inertia can be subtracted from the total, leaving only the rotational inertia of the ring.
2. Take the ring off the apparatus and repeat the steps under "Finding the Acceleration of the Ring and Disk" for the disk alone. Record the results in Data Table 3.

Calculations

3. Calculate the experimental value of the rotational inertia of the ring and disk together, and record the value in Data Table 3.

Wireless Rotary Motion Sensor Experiment 2: Rotational Inertia of Disk and

4. Calculate the experimental value of the rotational inertia of the disk alone and record the value in Data Table 3.
5. Subtract the rotational inertia of the disk from the total rotational inertia of the ring and disk, and record this as the rotational inertia of the ring alone.
6. Use a percent difference to compare the experimental values to the theoretical values.

Data Table 3: Results

Item	Rotational Inertia
Ring and Disk	
Disk alone	
Ring alone	
Percent difference for the disk	
Percent difference for the ring	

Experiment 3: Conservation of Angular Momentum

Equipment Required*	Equipment Required*
Wireless Rotary Motion Sensor (PS-3220)	Base and Support Rod (ME-9355)
PASCO Data Collection Software	Mass and Hanger Set (ME-8979)1
Mini-Rotational Accessory (CI-6691)	Triple Beam Balance (SE-8723)
Calipers (SF-8711)	Paper clips (for masses <1g)

*See the PASCO Web site at www.pasco.com for more information

Purpose

A non-rotating thin ring is dropped onto a rotating disk, and the final angular speed of the system is compared with the value predicted using the principle of the conservation of angular momentum.

Theory

When the ring is dropped onto the rotating disk, there is no net torque on the system since the torque on the ring is equal and opposite to the torque on the disk. Therefore, there is no change in angular momentum; angular momentum (L) is conserved.

$$L = I_i \omega_i = I_f \omega_f$$

where I_i is the initial rotational inertia and ω_i is the initial angular speed of the disk and I_f is the final rotational inertia and ω_f is the final angular speed of the disk and the ring together.

The rotational inertia of a disk is given as:

$$I_i = \frac{1}{2} M_1 R^2$$

and the final rotational inertia of a disk and ring together is:

$$I_f = \frac{1}{2} M_1 R^2 + \frac{1}{2} M_2 (r_1^2 + r_2^2)$$

where M_1 is the mass of the disk, M_2 is the mass of the ring, R is the radius of the disk, and r_1 and r_2 are the inner and outer radii of the ring.

Based on this, the final rotational speed is given by:

$$\omega_f = \frac{M_1 R^2}{M_1 R^2 + M_2 (r_1^2 + r_2^2)} \omega_i$$

Setup

1. Mount the Wireless Rotary Motion Sensor to a support rod. Remove the thumbscrew from the 3-Step Pulley. Place the disk directly on the square-shaped top of the 3-step Pulley as shown in Figure 3.1.
2. Place one of the disk alignment guides on the disk. Replace the thumbscrew to the 3-Step Pulley so that the disk alignment guide is firmly held in place.

Wireless Rotary Motion Sensor Experiment 3: Conservation of Angular Momen-

3. In the data collection software, create an experiment to measure the angular velocity (in radians per second) versus time (in second) of the disk before and after the ring is dropped on top of it.
 - In PASCO Capstone, for example, drag the Graph icon from the Displays palette to the workbook. Select "Angular Velocity (rad/s)" for the vertical axis, and "Time (s)" for the horizontal axis.

Procedure

1. Hold the thin ring just above the disk alignment guide that is on the top of the disk.
2. Give the disk a spin with your hand and click Record to begin recording data.
3. After about 25 data points have been recorded, drop the thin ring onto the spinning disk. See Figure 3.1.
4. Click Stop to end data recording after the disk and ring have made a few rotations.
5. In the Graph display, select the region of the data that represents when the ring was dropped onto the disk.
6. In the Graph display, select the data analysis tool that shows the coordinates of any point in the plot of data and move the cursor to the data point that is immediately before the collision. Record the Angular Velocity at this point as the initial angular velocity in the Data Table.
7. Move the cursor to the data point immediately after the collision. Record the Angular Velocity at this point as the final angular velocity in the Data Table.
8. Weigh the ring and disk and record their masses. Measure the inner and outer radii of the ring, and the radius of the disk. Record these values in the Data Table.


Figure 3.1 Setup for Conservation of Angular Momentum

Analysis

1. Calculate the theoretical value for the final angular velocity and record this value in the Data Table.
2. Calculate the percent difference between the experimental and theoretical values of the final angular velocity and record it in the Data Table.

Questions

1. How does the experimental result for the final angular velocity compare with the theoretical value for the final angular velocity?
2. What percentage of the rotational kinetic energy was "lost" during the collision? Calculate the energy lost and record the results in the Data Table.

$$\% \text{ KE lost} = \left(\frac{\frac{1}{2}I_i\omega_i^2 + \frac{1}{2}I_f\omega_f^2}{\frac{1}{2}I_i\omega_i^2} \right)$$


Data Table: Data and Results





Initial angular velocity	
Final angular velocity (experimental value)	
Mass of disk (M1)	
Mass of ring (M2)	
Inner radius of ring (r1)	
Outer radius of ring (r2)	
Radius of disk (R)	
Final angular velocity (theoretical value)	
Percent difference between experimental and theoretical values	
Percent of kinetic energy lost	

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